
Element

infinite
biomedical technologies

- | | | | |
|----|--|---|---|
| EN | | Installation and Setup Guide for the Prosthetist | ➔ |
| DA | | Installations- og opsætningsvejledning til proteselæge | ➔ |
| DE | | Installations - und Einrichtungsanleitung für den Orthopädietechniker | ➔ |
| ES | | Guía de instalación y configuración para el protésico | ➔ |
| FI | | Proteesintekijän asennus- ja käyttöönotto-opas | ➔ |
| FR | | Guide d'installation et de configuration pour le prothésiste | ➔ |
| NL | | Installatie- en configuratiehandleiding voor de prothesemaker | ➔ |
| NO | | Installasjonsveiledning | ➔ |
| PT | | Guia de Montagem e Ajuste para o Ortoprotésico | ➔ |
| SE | | Installations- och konfigurationsguide för protesläkaren | ➔ |

Element

Installation and Setup Guide
for the Prosthetist

infinite

biomedical technologies

All rights reserved. Element is a trademark of Infinite Biomedical Technologies, LLC.

This document provides information for the prosthetist who will be installing Element and IBT Electrodes.

Contains FCC ID: XDULE40-D2


Contains IC: 8456A-LE4D2

Infinite Biomedical Technologies, LLC.
8 Market Place, Suite 500
Baltimore, MD 21202
(443) 451-7175
www.i-biomed.com
SRN: US-MF-000007619

www.i-biomed.com/support.html

Mdi Europa GmbH
Langenhagener Str. 71
30855 Langenhagen
Germany
SRN: DE-AR-000006218

LEGEND OF SYMBOLS USED

Medical Device

Consult Instructions for Use

Single Patient, Multiple Use

Keep Dry

Distributor

Importer

CE Mark

Serial Number

Translated

European Authorized Representative

Manufacturer

Catalogue Number

Model Number

This symbol is used throughout the guide to indicate important cautionary information. Text following this symbol should be read carefully.

R_x Only

Caution: Federal law restricts this device to sale by or on the order of a prosthetist.

This device includes an RF transmitter or applies radio frequency electromagnetic energy.

Element

INSTALLATION AND SETUP GUIDE FOR THE PROSTHETIST

Thank you for choosing Element to provide EMG inputs for an upper limb myoelectric prosthesis. This guide will familiarize you with Element and help you install the hardware and software.

Have any questions? We're happy to help. Call us or send us an email.

(443) 451-7175

support@i-biomed.com

INTENDED USE

The Element system is intended to detect, process, and transmit physiological signals for use with a prosthesis

INDICATIONS FOR USE

The Element system is intended to be used exclusively for myoelectric exo-prosthetic fittings of the upper limb.

CONDITIONS FOR USE / PATIENT TARGET GROUP


Element is intended for use on one patient only, for users with unilateral or bilateral amputation, hand, forearm and upper arm amputation or dysmelia.

Use of the product by another person is not approved by the manufacturer. Installation of the system should be performed exclusively by a licensed prosthetist or technician. Any unauthorized handling or installation of Element could void its warranty.

Table of Contents

1 Meet Element _____	5
2 Component Description _____	7
3 Specifications _____	8
4 Installation _____	9
Before you Begin _____	9
Connecting IBT Electrodes and Batteries _____	10
Powering Element with FlexCell _____	11
Element Software _____	12
Incorporating IBT Electrodes Into Socket _____	17
Installing Element into the Prosthesis _____	21
5 Testing Element _____	26
Troubleshooting _____	26
6 Maintaining Element _____	26
Maintenance _____	26
Disposal _____	27
Repairs, Returns, and Warranty _____	27
7 Safety and Warnings _____	28
8 Regulatory Info _____	29

1 Meet Element

The Element System has three components:

- IBT Electrodes (up to 2)
- Signal Processing Box
- Element Desktop Software

Element should only be powered by FlexCell Batteries.

Element outputs standard envelope EMG signals detected from the IBT Electrodes placed on the user's skin. The Element System is an alternative to standard suction socket myoelectrodes, with the advantages of a lower profile, digital signal processing, and wireless gain adjustment. Element does not replace or modify any functionality of connected prosthetic components.

Element is compatible with most hands, wrists, and elbows that accept standard analog EMG electrode inputs. The Element system is typically sold with three-port kidney-style connectors to connect with quick-disconnect terminal devices. Alternative connectors may be available upon request. If you have a question about compatibility with other devices, please contact us.

For detailed information on connections and cables, refer to the Element Connections Guide.

The Element system was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks to the wrist unit (pushups, downhill mountain biking) or extreme sports (free climbing, paragliding, etc.). Furthermore the Element system should not be used for the operation of motor vehicles, heavy equipment (e.g. construction machines), industrial machines or motor-driven equipment.

2 Component Description

IBT ELECTRODES

The electrodes detect and amplify raw electromyography (EMG) signals from the user's skin. The electrodes plug into the signal processing box.

SIGNAL PROCESSING BOX

The signal processing box collects and filters the electrode EMG signals, and outputs envelope EMG signals to the terminal device. The signal processing box contains a Bluetooth module, which allows Element to communicate with the desktop software.

ELEMENT SOFTWARE

The user can visualize EMG signals and adjust electrode gains through the desktop software.

3 Specifications

Dimensions (Element Box LxWxH)	38mm x 22.8mm x 3.85mm
Dimensions (IBT Electrodes LxWxH)	28.8mm x 16.8mm x 6.7mm
Temperature range (use)	-10°C to +50°C (14°F to 122°F)
Temperature range (transport/storage)	-20°C to +65°C (-4°F to 149°F)
Humidity range (use)	45% - 75%
Humidity range (transport/storage)	15% - 93%
Atmospheric pressure range	860 hPa - 1060 hPa
Input voltage	5 to 10V
Maximum Output Current	3A
Compatible battery	FlexCell
Recommended battery capacity	Depends on terminal device. Contact us for recommendations.
Expected service life	3 years
Compatible electrode	IBT Electrodes
Bluetooth	FCC, IC, CE, RoHS, and Bluetooth® 4.0 Certified ISM 2.4GHz module

For FlexCell

Voltage Output	7.4V DC
Capacity Range *	550 mAh - 2200 mAh
Current Output	Up to 7A
Temperature range (use)	0°C to +49°C (32°F to 120°F)
Temperature range (transport and storage)	0°C to +49°C (32°F to 120°F)

* Capacity range is dependent on how many FlexCell batteries are installed.

4 Installation

BEFORE YOU BEGIN

Included in the Package

- Element signal processing box
- IBT Electrodes
- Molding dummies for IBT Electrodes
- Molding dummy for signal processing box (if requested)
- USB Thumb drive containing Element desktop software
- Bluetooth Adapter
- FlexCell batteries (if ordered with Element)

What You'll Need

- PC
- FlexCell batteries (if not ordered with Element)
- Terminal device
- Coaxial plug (if not using a wrist or elbow)
- Lamination collar parts


Terminal Devices That Have Been Tested For Compatibility With Element

Hands		Wrists	Elbows
SensorHand Speed	bebionic small	MC ProWrist (with 4 or 6 band coaxial plug)	DynamicArm
MyoHand VariPlus Speed	i-limb access		
ProHand	i-limb ultra		
ProETD	i-limb ultra revolution	Ottobock OB 10S17 with Myorotronic	
Steeper MyoHand	i-limb quantum		
bebionic3			

CONNECTING IBT ELECTRODES AND BATTERIES

WARNING: Turn the battery OFF before plugging in any cables. Do not attempt to connect the battery or electrodes in a way that is not specified. This could cause permanent damage to Element, the battery or electrodes.

POWERING ELEMENT WITH FLEXCELL

Turning FlexCell On and Off


To turn on FlexCell, press and hold the charging port button for 3 seconds. To turn off FlexCell, press and hold the Control Unit button for 3 seconds.

For detailed information on additional FlexCell functions, please refer to the FlexCell Instructions for Use.


ELEMENT SOFTWARE


1


3


2


4

5

6


Adjusting Gains

Click the "Gain Settings" button to hide or show the gain sliders. To change the gains, drag the slider button left or right to adjust the gain for each electrode.


If the gain is modified, the new gain is sent automatically to the signal processing box and the **changes are immediately reflected** in the signal viewer. However, to **save these gains permanently, you must press the "Save Gains" button**, which becomes red after any gain is adjusted. If you try to close the application without saving the gains, a popup window will ask if you want to save the gains or discard the changes.

The figure above shows signals in single plot mode. The gains were changed but not saved, therefore the Save Gains button is active. The signals are currently active and scrolling.


The figure above shows the signals in multi-plot mode. Gains were changed and saved, therefore the Save Gains button is inactive. The signals are also paused and the Play Graph button can be clicked to reactivate them.

INCORPORATING IBT ELECTRODES INTO SOCKET

Selecting Ideal Electrode Sites

Site selection for the IBT Electrodes should be performed using the same techniques for standard dry electrodes. Electrode signals can be viewed through the Element software to determine contraction strength in various locations.

NOT TO SCALE

WARNING: Ensure IBT Electrodes are placed only on intact skin.


WARNING: After removing the prosthesis, always check for skin irritation and pressure sores under the electrode sites. If irritation is present, discontinue the use of the prosthesis until the skin is healed and check if socket adjustments are necessary to alleviate pressure.


Fabricating Electrodes into the Inner Socket

Once electrode sites have been marked on the plaster cast, place the electrode dummies on the cast and check to make sure they lay flat against the plaster surface. Make sure the surface of the electrode dummy marked "CAST SIDE" is placed against the plaster cast (see image).

For more detailed instructions, refer to the Electrode Fabrication Tips guide.

Secure 2 of the nails provided (14ga or 15ga standard nails, roughly 1/2" long) through the two holes in the electrode dummy to the cast. Alternatively, you may use thin double stick tape to secure the dummy to the surface of the cast.


Inner Socket Fabrication


IBT Electrodes are designed to fit in a self-retaining pocket that is formed into the inner socket using the provided electrode molding dummies. The same dummies can be used for vacuum-forming thermoplastic inner sockets and rolled silicone inner sockets.

For more detailed instructions, refer to the Electrode Fabrication Tips guide.


Inner Socket Fabrication - Thermoplastic


Inner Socket Fabrication - Silicone

Installing the Electrode Into the Inner Socket


Note: If a suction socket is being created, completely sealing the electrode cable exit helps maintain consistent suction. The cable exit can be sealed with silicone adhesive. To seal the cable exit with silicone adhesive, apply a dime-sized amount of it to the back of the electrode and about 5mm up the cable protruding from the electrode. Do this before snapping the electrode into the pocket. Wait 15-20 minutes for this to dry before testing suction.

NOT TO SCALE

Outer Frame

INSTALLING ELEMENT INTO THE PROSTHESIS


Element has the following electrical output characteristics:

- EMG Signal (5V Envelope)
- Element is a pass-through for the battery, therefore output voltage matches input battery voltage (5-10V)

Element is compatible with all terminal devices that accept the aforementioned electrical outputs. Contact IBT for an up-to-date list of all compatible devices.


WARNING: For all setups below, ensure the battery is OFF before plugging in any cables. Do not attempt to connect the Element output cables in a way that is not specified. This could cause permanent damage to the hand and Element.

Hand Only

Hand and Wrist

Ottobock Wrist (10S17) with MyoRotronic (13E205)

CABLE HAS BLACK STRIPE

CABLE HAS NO BLACK STRIPE


CABLE HAS NO BLACK STRIPE

CABLE HAS BLACK STRIPE


WARNING: Make sure that the Element red connector is plugged into the coaxial plug power pins labeled "1". Incorrectly connecting the red connector could cause permanent damage to the hand and Element.

Motion Control ProWrist (P/N 5010056 - 5010058)

Elbow

Ottobock DynamicArm

5 Testing Element

Please test Element before fitting the patient with the prosthesis. Connect all parts together and power ON the prosthesis. Connect to the software and tap on each electrode. You should observe peaks on the signals and the terminal device should move.

Note: The patient's first use of the device must be in the presence of a trained practitioner.

TROUBLESHOOTING

For more detailed information, refer to the Element Troubleshooting guide.

6 Maintaining Element

MAINTENANCE

Cleaning Electrodes

For best signal responses from your EMG electrodes, we recommended cleaning the electrodes regularly to remove sweat residue and maintain hygiene. Clean the surface of each electrode and the inner socket with a mild soap solution as needed, at least a couple of times a week.

WARNING: DO NOT SUBMERGE whole electrode or arm in a bath of water. This could permanently damage the electrodes.

WARNING: DO NOT USE harsh chemicals, such as acetone, bleach, kitchen cleaners, etc. This could permanently damage the electrodes and the skin.

IBT Electrode and FlexCell Battery Replacements

IBT Electrodes and FlexCell batteries can be easily replaced if they are damaged or stop working. If you need to replace the battery or electrodes, simply unplug the battery or electrodes from the Element input connectors. Follow the instructions in Section 4 to plug in the new electrode or battery. For any abnormal issues, discontinue use and contact IBT for support.

DISPOSAL

The Element system components should not be thrown away with common household waste. Dispose of Element components by either returning the unit to IBT or taking the unit to an official electronics disposal site.

REPAIRS, RETURNS, AND WARRANTY

Please contact IBT at support@i-biomed.com regarding repairs and returns. The Element signal processing box and IBT Electrodes come with a 1-year manufacturer's defect warranty. Details of the warranty are enclosed separately.

7 Safety and Warnings

Liquid damage warning: Element is not waterproof. Please advise the user to avoid submerging or spilling liquid on or into their prosthesis. This may cause Element to be permanently damaged. If the user will be wearing their prosthesis in wet environments, ensure that Element is sufficiently protected from the external environment.

Battery use warning: Element must be powered by a FlexCell battery system from IBT. We have conducted safety testing to confirm compatibility between Element and FlexCell batteries.

Modification WARNING: Any unauthorized modification to Element can pose a safety risk to the user and will void the warranty. Changes or modifications not expressly approved by Infinite Biomedical Technologies, LLC could void the user's authority to operate the equipment.

WARNING: DO NOT USE harsh chemicals (such as bleach, kitchen cleaners, etc.) to clean electrodes. This could permanently damage the electrode. If you are not sure if a specific chemical is considered harsh, please call us prior to cleaning the electrode.

WARNING: Do not use Element during safety critical tasks.

WARNING: IBT Electrodes are only to be used with the Element signal processing box or other compatible IBT systems.

Element needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this document.

Please note that portable and mobile RF communications equipment can affect Element.

WARNING: Use of accessories and prosthetic parts other than those specified in the product order may result in malfunction of the Element system and can cause increased emissions and decreased immunity of the Element system.

INCIDENT REPORTING

Any serious incident needs to be reported to Infinite Biomedical Technologies and to the competent authority of the EU member state in which you and that patient are established.

IBT can be reached at:

Infinite Biomedical Technologies, LLC.
8 Market Place, Suite 500
Baltimore, MD 21202, USA
Phone: +1 (443) 451-7175
E-mail: info@i-biomed.com
Website: www.i-biomed.com

Contact information for competent authorities can be found at:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

8 Regulatory Info

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device. (Continued on next page).

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This product has been tested and verified to ensure that there are no issues or concerns regarding reciprocal interference. This includes EMI, EMC and RF.

This product has been certified and tested by 3rd party testing facilities to the following standards:

IEC 60601-1, 3rd Edition
IEC 60601-1-2, 3rd and 4th Edition
IEC 60601-1-11, 1st Edition
IEC 61000: See next page

Also compliant as per CISPR 11:2015

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The Element System is intended for use in the electromagnetic environment specified below. The customer or the user of the Element System should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 2	The Element System must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF Emissions CISPR 11	Class A	The Element System is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	


Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The Element System is intended for use in the electromagnetic environment specified below. The customer or the user of the Element System should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	Not applicable
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	Not applicable	Not applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the A.C. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The Element System is intended for use in the electromagnetic environment specified below. The customer or the user of the Element System should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the Element System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance Not applicable
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 Ghz	10 V/m	$d=0.35 \sqrt{P}$ 80 MHz to 800 MHz $d=0.7 \sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Element System is used exceeds the applicable RF compliance level above, the Element System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Element System.			

Recommended separation distances between portable and mobile RF communications equipment and the Element System			
The Element System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Element System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Element System as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{V_i} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{E_i} \right] \sqrt{P}$	800 MHz to 2,5 GHz $d = \left[\frac{7}{E_i} \right] \sqrt{P}$
0,01	Not applicable	0.035	0.07
0,1	Not applicable	0.1106	0.221
1	Not applicable	0.35	0.7
10	Not applicable	1.106	2.21
100	Not applicable	3.5	7
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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
Stepes
535 Mission Street, 15th Floor
San Francisco, CA 94105
USA

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